PECG HC3A250
User Manual

1.Introduction:
HC3A250 is a miniature wearable dynamic PECG recorder. HC3A250 may be used to record from 1 minute to 48 hours PECG (or Holter monitoring), HC3A250 also has alert message function with alert message triggered by either panic button push, or heart beat stop. The HC3A250 should not be used for newborn baby weight less than 10kg. HC3A250 does not display pulse, or diagnostic function like heart rate algorithm and ST segment recognition function. The professional analysis ECG report (additional charge) would be provided by a third party.

Product basic information:
1)Product Name: PECG
2)Model: HC3A250
3)Production date: see the backside of the PECG
4)Expiry date: 3 years
5)Power: built-in rechargeable polymer lithium battery
Nominal voltage: 3.7V
Capacity 400mAh (1.48Wh)
6)Products and accessories:
a)1pc X PECG
b)12pcs X Single-use ECG electrodes
c)1 x USB data cable (connecting from PC to PECG)
d)1 x OTG USB Cable (connecting from mobile to PECG)
e)1 x Warranty card / user manual
f)1 x Certificate

Structure Description:

2.How to use :
(1)Preparing for device with electrodes
Put PECG in power-off state, assembly sticky patches called electrodes and remove the electrode protective films

(2)Wear PECG onto the chest
Prior to application of electrodes, suggested to shave or clean the body hair at electrode area (red oval area), thoroughly clean oily skin on the chest where the electrodes are placed, and keep the operation key of the PECG upward, and apply electrodes firmly to the skin. Make sure all electrodes are placed at site specified in the figure below, two electrodes in lower part of PECG should be placed on both sides of the chest wall where is nearly epigastric site.

(3)Open PECG to start recording
Switch on PECG, press & hold button till green led light on, after a couple of seconds, the green light will flash and PECG data collection starts immediately.

(4)Close PECG to stop monitoring
Switch off PECG, press & hold button continuously till green light on, after a few seconds, a green light flash turns off automatically, PECG data collection stops (normal use continuously 24 hours).

(5)Remove PECG from chest
When data collection stops, press ECG with a hand (see below diagram), remove electrodes from the chest slowly from outside to inside, and then pull out PECG & dispose of used electrodes into garbage bin. DON'T force to tear away electrodes from your chest to avoid any skin damage.

(7)Emergency Call: Push the “Switch” twice, yellow LED light will flash, the buzzer 3 times consecutively in 2 seconds, at the same time a short message would be sent to the emergency contact(s) via APP service. When push the Switch again, the buzzer stops the call.

(8)Heart Beat Stop Alarm: If the user occurs cardiac arrest heart beating for few seconds, will enter into an alarm state, PECG will trigger an emergency call for help, yellow LED light will flash, the buzzer sound 3 times consecutively in 2 seconds, at the same time a short message would be sent to the emergency contact via APP service. When push the Switch again, the buzzer stops the call.

3.Function Description:
PECG status description
PECG has three states: working, shut-down and charging / USB connection. When the PECG is in the shut-down state, press the “Switch” for about 3 seconds and the green light will be on for about 5 seconds, then PECG will enter the working state. When the PECG is in working state, press the Switch key for about 3 seconds. The green light will be on for 3-5 seconds and the PECG will enter in shut-down state. When the PECG is in working state or shut-down state, connect the USB data cable to the computer, which can enter the charging and USB connection state. When HC3A250 is in the charging / USB connection state, ECG data can not be collected. The following table is a detailed description of the three states:

Status name	Status description	Indicator light, buzzer reminder
Working state	When the PECG electrode is in good contact with the skin, the green light flashes once in 2 seconds;	When the PECG electrode is in good contact with the skin, the green light flashes once in 2 seconds;
	When the PECG electrode is in bad contact with the skin, the yellow light flashes once for about 2 seconds;	When the PECG electrode is in bad contact with the skin, the yellow light flashes once for about 2 seconds;
	The buzzer will sound once in 2 seconds when the battery is low.	The buzzer will sound once in 2 seconds when the battery is low.
	Once user heart beat stopped for few seconds, the PECG will trigger an “Emergency Call”, the buzzer will sound 3 times in 2 seconds. When push the switch again to cancel the “Emergency Call” & the buzzer sound.	Once user heart beat stopped for few seconds, the PECG will trigger an “Emergency Call”, the buzzer will sound 3 times in 2 seconds. When push the switch again to cancel the “Emergency Call” & the buzzer sound.
Shut-down state	Entering this state, ECG data collection stops	—
Charging / USB connection state	Entering this state, indicates that the PECG starts charging or charging while transferring data to the computer	During charging, the yellow light flashes once in 2 seconds. The green light stays green when fully charged or charging is complete.

The three state switching modes are shown below:

```
graph TD
    Working[Working state] -- "Press and hold the key in 3 seconds" --> Shut-down[Shut-down state]
    Shut-down -- "Press and hold the key in 3 seconds" --> Working
    Shut-down -- "Connect the PECG to the PC/mobile phone with a USB data cable" --> Charging[Charging / USB connection state]
    Charging -- "Disconnect the PECG from the PC/mobile phone's USB connection" --> Shut-down
```

Cautions :
◆ Please read the instructions carefully before using the PECG.
◆ Keep the PECG away from heat, moisture & water.
◆ Electrode is for one-time use only. Repeated use may cause skin infection.
⊗ Please do not disassemble and replace the battery, if necessary, please contact our authorized after-sales service department or dealer for processing.

4.Q & A:

Problem	Causes	Suggestions
Noisy waveform	Electrode does not make good contact with the skin caused by sweat, misplacement or reuse electrode	Make sure the chest skin is dry, clean and replace with new electrode
The computer does not detect the flash drive on the PECG after connection.	The computer may not have USB driver	Install USB driver in PECG
PECG device cannot connect with computer	USB cable may be damaged or does not have a good contact with computer	Use a new USB cable, pull out the cable from computer and re-connect again.
"TICK" sound from PECG	The battery is at a low power	Recharge the PECG device

Symbol:

	Fail to observe the warning may cause seriously injury or death to the user. May also cause the PECG to partial or completely damaged .
	Fail to observe the note may cause injury to the user. May also cause the PECG to work improperly.
	CF application part

5.Disclaimer:

● All data collected by PECG can not be used for medical diagnosis unless recognized by BISA authorized medical units, experts or doctors.
● The company will not take any responsibility for the misdiagnosis, misleading and loss caused by the self interpretation of the ECG data.
● The risk assessment based on ECG signals collected by ECG is not perfect. We will not take any responsibility for any loss caused by the Risk assessment which may affect user's health condition.
● The Company expressly disclaims any liability for loss of any direct, indirect or consequential damages caused by the PECG, including but not limited to loss of revenue, loss of profits, business interruption or loss of data, regardless of any form of action or legal claim of the liability, and even including any prior notice of any loss incurred.

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